




CNT240 Digital Pathology Slide Scanner

Digital Pathology: Accuracy in a Single Scan
All-in-One Solution

- Fully Automated
- Multifunctional
- High-Throughput



CNT240

Digital Pathology Slide Scanner


Integrates multiple functions to create a standardized slide scanning workstation; a one-stop solution for digital pathology



CORE STRENGTHS



STAINING EFFECT



TECHNICAL PARAMETER

Model	CNT240
Throughput	240 pieces
Scanning Speed	No less than 4 cm ² /min
Image Resolution	0.25um/pixel(±15%)
Objective Lens	20x, 40x real-object objectives
Scanning Camera	12.3 megapixels
Preview Camera	2 megapixels
Brightfield Illumination	High-power LED
Image Format	SVS / JPG / HDS / TIFF / BMP
Software	Automated Digital Microscope Scanning System